

SAMYANG

# croquis

SYB-MD-CRO-202411V1

SAMYANG

295 Pangyo-ro, Bundang-gu, Seongnam-si Gyeonggi-do 13488, Korea  
croquis@samyang.com | www.syb-smart.com



# croquis®

## Samyang Holdings Corporation

Business Areas .....	4
Major Achievements	
Major Location .....	5

## croquis®

Brand Motive .....	8
MONO .....	10
SCREW .....	12
VOLUME .....	14
BARBED CUTTING .....	16
BARBED MOLDING .....	18

## SYB Aesthetics

Art Technology .....	6
----------------------	---

## Guaranteed Quality

Storage Test .....	20
BHS Test .....	21
* BHS (Barb Holding Strength)	

## Brand Vision

SMART Academy .....	22
(Medical Aesthetics Education Platform)	




# Samyang Holdings Corporation

Through accumulated experience and technology

With more than 25 years of experience in suture manufacturing and polymer technology, Samyang Holdings now intends to expand its portfolio for Medical Aesthetics.

## Business Areas



## Major Achievements



## Major Location



Samyang Discovery Center  
(Headquarter & R&D Center)

Completion	2016
Related Department	Marketing & Sales / Regulatory Affairs / Planning / R&D Center
Location	Seongnam-si, Gyeonggi-do, Korea

### Pharmaceuticals Plant



**COMPLETION**  
1997

**PRODUCT**  
Injection, API

**LOCATION**  
Daejeon, Korea

**CERTIFICATIONS**  
EU GMP, PMDA (Japan), GMP, KGMP (Korea)

### Medical Devices Plant




**COMPLETION**  
1996

**PRODUCT**  
Sutures, Hemostats, Patches

**LOCATION**  
Daejeon, Korea

**CERTIFICATIONS**  
KGMP, CE Mark, ISO13485 (Medical Devices), KGMP (Pharmaceuticals)

### Samyang Biopharm USA Inc.




**COMPLETION**  
2018

**BUSINESS MODEL**

- NRDO (No Research Development Only)
- Therapeutics for Oncology and Rare Diseases
- License-in / Out

**LOCATION**  
Cambridge, MA, USA




## SYB Aesthetics

Everyone is born beautiful.

The definition of beauty has changed from wanting to look like someone to bringing out your beauty within. Samyang believes in value as art. SYB Aesthetics is derived from combining art with technology. By bringing out your embedded beauty with Samyang's technology, SYB Aesthetics provides a helpful and safe solution to stay young and beautiful.


### Art Technology





# croquis®

## Brand Motive

croquis® is motivated by the drawing technique of describing an object in a short period of time. Find your embedded beauty through skin rejuvenation and tissue repositioning.



# RETURNING TO THAT YOUTH



Tightening ●●● | Volume ●○○○

croquis<sup>®</sup> MONO stimulates neocollagenesis, which not only improves your overall skin quality but also fully absorbs into the body over time.

Material : PDO (Polydioxanone)  
 Contains : 60ea / box (10ea / Al pouch \* 6)  
 Thread Type : \_\_\_\_\_

Needle		Thread	
Gauge	Length (mm)	USP	Length (mm)
25G	60	5 - 0	90
	90		150
27G	60	6 - 0	90
	40		50
29G	50		70
	25		30
30G	25		30
	40		50



Tightening ●●● | Volume ●●○

croquis<sup>®</sup> SCREW stimulates neocollagenesis and improves your overall skin and wrinkle through volumization effect. SCREW type also gets absorbed completely into the body over time.



Material : PDO (Polydioxanone)

Contains : 60ea / box (10ea / Al pouch \* 6)

Thread Type :

Needle		Thread	
Gauge	Length (mm)	USP	Length (mm)
26G	50	5 - 0	70
	60		90
	90		130
27G	40	6 - 0	52
	50		70
29G	40	6 - 0	52
	50		70

# B RINGING THAT SMOOTH AND TIGHT SKIN



 VOLUME

Tightening ●●● | Volume ●●●

croquis<sup>®</sup> VOLUME effectively adds volume to deep wrinkles or fold, providing a smoother and more rejuvenated appearance.



Material : PDO (Polydioxanone)

Contains : 30ea / box (6ea / AI pouch \* 5)

Thread Type : 



L-Cannula		Thread	
Gauge	Length (mm)	USP	Length (mm)
19G	40	6 - 0 x 9ea	60
	60		90
21G	40	7 - 0 x 9ea	60
	60		90

# A DDDING THAT VOLUME TO Y OUR SKIN





# PULLING AND HOLDING THAT BEAUTY



## BARBED CUTTING

Skin Quality ●●●

Lifting ●●●

croquis<sup>®</sup> BARBED CUTTING is used for tissue fixation and repositioning. With absorbable component and neocollagenesis, BARBED CUTTING type provides a younger and beautiful look.

Material : PDO (Polydioxanone)

Contains : 30ea / box (6ea / Al pouch \* 5)

Thread Type :   
 (3D, Bi-Directional) (Reverse Bi-Directional)

L-Cannula		Thread		Barb
Gauge	Length (mm)	USP	Length (mm)	
18G	100	2	150	(3D, Bi-Directional)
19G	60	0	105	
	90		140	
20G	60	2 - 0	105	(Reverse Bi-Directional)
	90		140	
19G	40	2	50	
	50		70	

**BARBED MOLDING**

Skin Quality ●●● | Lifting ●●●+

croquis<sup>®</sup> BARBED MOLDING is used for a stronger lift. By using the pine tree-like barb shape, BARBED MOLDING type gives a helpful solution to skin laxity.


Material : PDO (Polydioxanone)

Contains : 30ea / box (6ea / Al pouch \* 5)

Thread Type :  (2D, Bi-Directional) (2D, Multi-Directional)

L-Cannula		Thread		Barb
Gauge	Length (mm)	USP	Length (mm)	
18G	100	2	150	● (2D, Bi-Directional)
19G	100	1	185	
19G	70	1	155	● FIX (2D, Multi-Directional)
	100		185	
20G	70	0	155	
	100		185	

# EN HANCING AND KEEPING THAT LIFT



# Guaranteed Quality

High stability of croquis<sup>®</sup> can solve quality-related issues during storage and distribution processes. Such stability provides mutual trust to global customers.

## 1) Storage Test

Result after 5-year accelerated aging test with competing brands


**MONO**  
SCREW  
TYPE : MONO  
USP : 6 - 0



**BARBED CUTTING**  
USP : 3 - 0, 4 - 0



**BARBED MOLDING**  
USP : 2 - 0




## 2) BHS Test

croquis<sup>®</sup> BARBED CUTTING type's high BHS demonstrates effective facial lifting.

**BARBED CUTTING**


BHS : Barb Holding Strength (kgf)  
USP : 0



# Brand Vision

## croquis® as Market Standard

SYB Aesthetics provides high and safe quality products through well-controlled and strict manufacturing process.



# Education Platform

## SMART Academy

SMART Academy is an online training platform for medical aesthetics experts.

- Product Information
- Guide Protocol
- e-Course: Procedural Video



Samyang Medical Aesthetics Research & Training

Draw your beauty,  
Capture the moment

croquis®



croquis<sup>®</sup>

samyang



syb-aesthetics



SMART Academy



E-brochure

